



**SCW MEDICATH LTD.**

No.4 Baolong 6<sup>th</sup> Road, Baolong Industrial Town,  
Longgang District, Shenzhen, 518116, Guangdong, China  
TEL: 86-755-89312160 FAX: 86-755-89312239

## Product Summary Sheet

Reference: 43.01.20004

Description: Postpartum Balloon 24F, 55cm, 500mL, with Rapid Instillation

### Content


Component		Qty/Unit
Postpartum Balloon	Release channel	1
	Catheter	1
	Balloon	2
	Tip	1
	Screwed joint	1
	Two-way cock	1
150cm connecting tubing with puncture device	Connecting tubing	1
	Puncture device	
Syringe (60mL)		1

### Photo



### Packaging

Attribute	Material / Detail
Package	Tyvek® paper and Blister box
Quantity / Case	40 (4 Inner boxes of 10 )
Approx. Weight / Case	8.5kg
Shipper Dimensions	641*448*445mm
Sterilant	ETO



Peel To Open

SCW

Medicath  
Ltd.

Postpartum Balloon

(with Rapid Instillation Components)

550MM

Model:PPB-24F-II

Contents

postpartum balloon  
150cm extension tube with puncture device  
three-way cock with one-way valve  
60ml plastic injector

REF

LOT

SCW MEDICATH LTD.

ADD: No. 4, Baolong 6th Road, Baolong Industrial Town,  
Longgang District, Shenzhen 518116 Guangdong China  
TEL: (86)755-89312160/89312258  
FAX: (86)755-89312239  
EMAIL: sales@scwmed.com

EC REP

OBELIS S.A.

ADD: Bd. Général Wahis 53, 1030 Brussels, Belgium

CE

0197

SCW-LB-34

Available Version Number:A.3