

Product Specification Sheet

Glycerine Swabsticks 1115 & 1116



optimummedical®

● General information

General description: Pre-moistened swabsticks with water, glycerine, citric acid, sodium benzoate and lemon or blackcurrent flavouring.

Intended use: Pre-moistened swabsticks for oral refreshment.

Lifetime of device: Glycerine Swabsticks are intended to be discarded after single use. Not for long-term use.



Product code	Format
1115	Glycerine Swabsticks – lemon flavour
1116	Glycerine Swabsticks – lemon flavour

● Physical properties:

- Cotton swab diameter: 10mm ± 1mm
- Volume of liquid absorbed into cotton swab: 4.5ml ± 0.4ml
- Sterilisation: non-sterile
- Latex free
- Phthalate free
- Sugar free
- No materials of animal origin

● Packaging configuration:

3 swabsticks per sachet	65 x 145 mm [sachet]
25 sachets per inner box	160 x 85 x 145 mm [inner]
20 inners per outer box	440 x 350 x 310 mm [outer] Weight 7.5 Kg

All packaging is recyclable.

● Storage Conditions

- Store away from direct sunlight
- Keep dry
- Store at ambient temperature
- Shelf Life: 5 years
- Do not use if packaging is damaged

● Contains:

- Swabsticks: Glycerine, Water, Citric acid, Sodium benzoate, and Lemon or Blackcurrent flavouring absorbed into cotton wool.
- Polypropylene plastic handle.

