

Ultrasound Gel

The primary role of an ultrasound gel is to act as a conductor for the ultrasound waves, allowing them to pass into the body without being affected by the air. The secondary role of ultrasound gel during ultrasound procedures, is to allow the probe or medical device used for the procedure, to move smoothly, reducing resistance and potential discomfort to the patient.

IN THIS RANGE

Vue Ultrasound Gel	22
Sterile, single-use sachets	23
Non-sterile ultrasound gel	24

vue.TM

ultrasound gel

Vue Ultrasound Gel is a medical device designed for use during medical ultrasound procedures.

When selecting an ultrasound gel, it's important to select the right format of gel for the procedure and the patient.

- Classified as a Medical Device.
- Developed specifically for use during ultrasound procedures to have the ideal viscosity.
- Suitable for external ultrasound procedures.
- Sterile sachets suitable for use on high-risk patients.
- Water-soluble and non-sticky formulation.
- Won't stain or damage skin, clothing or medical equipment.
- Latex free.
- Available in sterile and non-sterile formats
- Available in a range of sizes and formats
- Available clear or blue.

Sterile, single-use sachets



20ml

Our single-use, sterile ultrasound gel formats are ideal for use in any healthcare setting where infection prevention and cross-contamination are key considerations.

Vue Ultrasound Gel double-wrapped sachets are recommended for use during procedures that require an aseptic technique.



Vue sterile, single-use sachets

PR CODE	PR NAME	SIZE	INNER PACK	OUTER PACK
1156S	Vue Ultrasound Gel	20ml	30	240
1156	Vue Ultrasound Gel	20ml double wrapped	30	240

Non-sterile ultrasound gel



250ml
1l

250ml and 1l bottles are ideal for use in a busy healthcare setting. Non-sterile ultrasound gel is suitable for use during low risk, general examinations, on intact skin.




Vue non-sterile ultrasound gel

PR CODE	PR NAME	SIZE	VARIATION	OUTER PACK
1157	Vue Ultrasound Gel	250ml bottle	Clear	40
1157B	Vue Ultrasound Gel	250ml bottle	Blue	40
1159	Vue Ultrasound Gel	1l bottle	Clear	12
1159B	Vue Ultrasound Gel	1l bottle	Blue	12

For more information on any
product in our range,
email enquiries@optimummedical.co.uk
or call +44 (0) 845 643 5479.

Our team are waiting to answer all your
questions and would be happy to arrange
for a free sample to be sent.

We also have a comprehensive range of catheterisation and catheter management medical devices?



For information on our **Ugo Urology** range
visit www.optimummedical.co.uk

For additional information
email enquiries@optimummedical.co.uk
or call +44 (0) 845 643 5479