

Catheterisation Jelly

It is recommended that the sterile lubricating jelly used during catheter insertion should be supplied in a pre-filled syringe, to allow effective instillation into the urethra, ensuring the entire urethra is lubricated before insertion of the catheter.

'The lubricant dilates and lubricates the urethra. It should be administered in a syringe to facilitate insertion of the gel into the urethra in accordance with local practice.'

(ENSH ² 2021)

IN THIS RANGE

Our catheterisation jelly range	14
OptiLube Active	15
OptiLube Active CHG Free	16
OptiLube	17

Our catheterisation jelly range

We believe in choice. Which is why we offer a range of catheterisation jelly pre-filled syringes, designed specifically for use during catheterisation procedures.

All our catheterisation jelly syringes are supplied sterile inside and out and designed to fall from their packaging onto a sterile field, supporting an aseptic technique.

- Each available in 6ml (female) and 11ml (male).
- Manufacturing date, expiry date and LOT number are printed on the tear off strips.
- Provided in a familiar and recommended syringe format.
- Easy-to-prime syringe for controlled application.
- Rounded tip for comfort when inserted into the urethra.

OptiLube Active



6ml
11ml

Containing local anaesthetic (Lidocaine 2%) and antiseptic (Chlorhexidine Gluconate 0.05%).



OptiLube Active

PR CODE	PR NAME	SIZE	INNER PACK	OUTER PACK
1160	OptiLube Active	6ml	10	200
1161	OptiLube Active	11ml	10	200

OptiLube Active CHG Free



6ml
11ml

Containing local anaesthetic (Lidocaine 2%). Ideally suited for patients with a contraindication or sensitivity to Chlorhexidine Gluconate (CHG).



OptiLube Active CHG Free

PR CODE	PR NAME	SIZE	INNER PACK	OUTER PACK
1167	OptiLube Active CHG Free	6ml	10	200
1168	OptiLube Active CHG Free	11ml	10	200

OptiLube Syringes



6ml

11ml

Containing no active ingredients. Ideally suited for patients with a contraindication or sensitivity to Lidocaine or Chlorhexidine Gluconate (CHG).



OptiLube

PR CODE	PR NAME	SIZE	INNER PACK	OUTER PACK
1125	OptiLube Syringe	6ml	25	250
1126	OptiLube Syringe	11ml	25	250