



Hydrophilic Coated Catheters

aloeoriginal

Our **Aloe Original** ISC catheter features our hydrophilic coating. When activated using water, the water-based coating ensures the entire catheter is lubricated, reducing friction during insertion and removal.

Aloe Original

Catalogue No.	Type	Colour	Charriere size	Diameter	Length	Pack size
ABP06	Paediatric	Green	6CH	2.0mm	28cm	30
ABP08	Paediatric	Blue	8CH	2.6mm	28cm	30
ABP10	Paediatric	Black	10CH	3.3mm	28cm	30
ABF08	Female	Blue	8CH	2.6mm	18cm	30
ABF10	Female	Black	10CH	3.3mm	18cm	30
ABF12	Female	White	12CH	4.0mm	18cm	30
ABF14	Female	Green	14CH	4.6mm	18cm	30
ABM08	Male	Blue	8CH	2.6mm	40cm	30
ABM10	Male	Black	10CH	3.3mm	40cm	30
ABM12	Male	White	12CH	4.0mm	40cm	30
ABM14	Male	Green	14CH	4.6mm	40cm	30
ABM16	Male	Orange	16CH	5.3mm	40cm	30
ABM18	Male	Red	18CH	6.0mm	40cm	30



aloeoriginal

KEY FEATURES:

- Hydrophilic coating activated using just water.
- DEHP Free.
- Nelaton tip.
- Coloured funnel for easy size recognition.
- Two smooth, rounded drainage holes for effective urine drainage and reduced discomfort during catheter insertion and removal.
- Packaging features a sticky tab to assist with positioning of the catheter during procedure set up.
- Suitable for draining directly into a toilet or receptable.
- Can be attached to a drainage bag if required.
- Available in paediatric, adult male and adult female sizes.

