

Non-coated Catheters


smoothnelaton

Our **Smooth Nelaton** is a non-coated, DEHP free ISC catheter with a Nelaton tip.

The **Smooth Nelaton** is an ISC catheter, which features smooth, rounded drainage holes to improve comfort during catheter insertion and removal.

Smooth Nelaton

Catalogue No.	Type	Colour	Charriere size	Diameter	Length	Pack size
SNP06	Paediatric	Green	6CH	2.0mm	28cm	30
SNP08	Paediatric	Blue	8CH	2.6mm	28cm	30
SNP10	Paediatric	Black	10CH	3.3mm	28cm	30
SNF08	Female	Blue	8CH	2.6mm	18cm	30
SNF10	Female	Black	10CH	3.3mm	18cm	30
SNF12	Female	White	12CH	4.0mm	18cm	30
SNF14	Female	Green	14CH	4.6mm	18cm	30
SNF16	Female	Orange	16CH	5.3mm	18cm	30
SNF18	Female	Red	18CH	6.0mm	18cm	30
SNM08	Male	Blue	8CH	2.6mm	40cm	30
SNM10	Male	Black	10CH	3.3mm	40cm	30
SNM12	Male	White	12CH	4.0mm	40cm	30
SNM14	Male	Green	14CH	4.6mm	40cm	30
SNM16	Male	Orange	16CH	5.3mm	40cm	30
SNM18	Male	Red	18CH	6.0mm	40cm	30



smoothnelaton

KEY FEATURES:

- Non-coated.
- DEHP Free.
- Can be lubricated using a water-based lubricant if required.
- Nelaton tip.
- Coloured funnel for easy size recognition.
- Two smooth, rounded drainage holes for effective urine drainage and reduced discomfort during catheter insertion and removal.
- Suitable for draining directly into a toilet or receptable.
- Can be attached to a drainage bag if required.
- Available in paediatric, adult male and adult female sizes.

