



Developed with clinician input, the **Ugo Pro** provides everything needed to support best practice insertion and removal of transurethral or suprapubic catheters in a community setting. Each pack contains a pre-wrapped sterile field and the essential routine components required to carry out the procedure, helping to minimise the risk of infection.

The **Ugo Pro** helps to minimise the unnecessary waste which would accumulate if each item were to be obtained separately and assists in streamlining stock management.

UGO PRO CONTAINS:

- 1 Clear instructions for use.
- 2 1 x outer wrapper – creates an aseptic field for the components to be placed upon.
- 3 1 x disposable apron.
- 4 1 x waste bag.
- 5 1 x grey waterproof drape – protects the bed from spilt fluids during procedure. Assists with aseptic technique.
- 6 1 x blue waterproof "O" drape - provides dignity during the procedure. Assists with aseptic technique.
- 7 1 x Ugo 10ml syringe – to remove fluid from Foley catheter balloon prior to catheter removal.
- 8 6 x non-woven swabs – to assist with removal and general catheterisation procedure technique.
- 9 5 x non-woven balls – for cleansing catheter insertion site during procedure (once coated in saline).
- 10 1 x 20ml sterile saline (0.9% w/v sodium chloride) - to empty over non-woven balls for cleansing catheter insertion site.
- 11 1 x clear tray with 2 compartments.
- 12 3 x pairs of medium nitrile **Ugo Pro** gloves – assists with best practice technique.

Ugo Pro (catheterisation procedure pack)

PR CODE	PR NAME	OUTER CASE
3400	Ugo Pro (Catheterisation procedure pack)	Case of 80

