



Easy to use, discreet and comfortable to wear, and designed in detail for safety and security. The **Ugo Catheter Valve** gives patients the confidence they need to carry on living independent lives. Using a catheter valve as an alternative to urine drainage bags may improve experience and outcomes for suitable patients.

The **Ugo Catheter Valve** is a small, sterile, single-use drainage valve which is used to allow intermittent urinary drainage from an indwelling urethral or suprapubic catheter.

The **Ugo Catheter Valve** features a universal stepped connector which is connected directly to the catheter, and a lever for controlled manual urine drainage. A **Ugo Catheter Valve** can remain in place for up to 7 days.

*"I use the **Ugo Catheter Valve** with my suprapubic catheter. It grips extremely well, and doesn't slip like some other valves which I have tried."*

**Christine P**, London, UK

### Ugo Catheter Valve

PR CODE	DESCRIPTION	INNER PACK	OUTER CASE
3000	Ugo Catheter Valve	Box of 5	Case of 200

## Key features:

- Offers a discreet alternative to urinary drainage bags for people assessed and deemed suitable.
- Allows the bladder to fill and hold urine.
- Can remain in place for up to 7 days.
- Suitable for connecting to a drainage bag for link drainage if required.
- Safety lever feature to prevent accidental opening.
- Universal stepped connector for secure fit
- Lever handle designed to allow easy single-handed operation.
- Contrasting colour safety lever for easy visual identification.

